# MDM万东

Precision & Vision



#### Beijing Wandong Medical Technology Co., Ltd.

Address: Building 3 No.9 Jiuxianqiaodong Road, Chaoyang District,100015,Beijing, P.R.China

Tel: +86-10-84575792/3/5/6

Fax: +86-10-84575794

Web: www.wandong.com.cn

E-mail: International@wandong.com.cn

Note: Design and specification subject to change without notice.



Panora Mew Oriental 1000LB

Full-frame Imaging without Stitching

# Full-frame Imaging without Stitching

The New Oriental 1000L series Panorama DR is equipped with an industry-leading 120cm flat-panel detector to create a high-definition and fast full-frame imaging platform without stitching, which brings a wide FOV to DR radiography, improves clinical efficiency, and adds another powerful tool for orthopedic diagnosis.

#### **Full View**

120cm length FPD
Full-frame imaging in one shot

#### MIX

Intelligent balanced imaging system

Consistent image quality

#### **SWA**

Super Wide-Angle Imaging System
Large FOV with high quality output

### **Multiple Clinic Applications**

Routine radiography/ Medical screening
Trauma First Aid / Orthopedics



01 02



# Full View 120cm Flat Panel Detector

Dedicated for clinical full-length imaging

Large field of view for rapid and accurate diagnosis.



Large FOV FPD






Monoblock









# No Stitching

Large FOV coverage, one exposure for whole limbs and vertebrae imaging.

## **Instant Imaging**

Instantaneous imaging of the whole limb or whole spine without complicated procedures, making the process simple and efficient.

#### **Low Dose Examination**


Full-limb, full-spine imaging in a single exposure, fast examination at low dose bringing more patient care.

03

# Extended Vision, Panorama Imaging

# **SWA Super Wide-angle Imaging System**

Restore the true image, getting what you see.



#### Flexible View

The field of view is automatically tracked and adjusted according to the anatomical programs.











120×43

# **MIX Intelligent Balanced Imaging System**

Combining multi-level filtration technology, multi-scale image decomposition algorithm and image fusion algorithm, breaking the body part limitation, allowing different body parts such as the skull, neck, chest, abdomen etc. to be imaged by one low dose shot, and in one uniform and consistent HD image.



















# Dedicated Care Everywhere


# **Applicable to Various Clinical Needs**

Explore multiple clinical applications and enable accurate clinical diagnosis.



# **Routine Radiography**

Daily radiography in various positions in the



# **Medical Screening**



#### Trauma First Aid

Wide range fast imaging Precise location of lesion





## 24hours online support

Quick response and remote technical support, on-line trouble shooting, spare parts inventory and quick delivery.

## Free software upgrade

Fully self-developed acquisition platform, free lifetime software upgrade, always providing you the excellent experience.

#### **Professional team**

Wandong imaging products has been installed in over 90 countries, more than 100 distributors and authorized service partners are ready to provide professional service.

# Regular maintenance

Regular maintenance by experienced and qualified engineers, guarantees the stable operation, prevents major failures, and extends the equipment life.



10